

Testing Equipments

CENTRIFUGES

For Routine Test


Model	C-854/6	
Max. Speed	RPM	3500
Max. RCF	'g'	1600
Rotor	S.O.H.	6 x 15ml

For Cross Matching


Models	R-303	R-304
Max. Speed	RPM	3500
Max. RCF	'g'	1250
Rotor	A.H.	8 x 15ml / 24 x 5ml

For General Purpose


Model	R-8 M	
Max. Speed	RPM	5250
Max. RCF	'g'	3600
Rotor: R-81 - 16 x 15 ml S.O.H.		

For Micro Hematocrite & Q.B.C.


Model	R-12M	
Capacity	ml	75mm x1mm
Tubes	No.	24 / 12
Max. Speed	RPM	18000
Max. RCF	'g'	18500
Rotor	: RM-1210	
Capillary	: RM-1212	
Reading Device	: RM-1213	

For RIA Test, Vacutainer tubes & Microtiter Trays


Model	PR-23	
Capacity	ml	108 x 5 ml
Tubes	No.	108
Max. Speed	RPM	7000
Max. RCF	'g'	7400
Rotor : R - 230 B with		
R - 233 B for RIA / Vacutainer		
R - 230 A for Microtiter tray		

SHAKER



For VDRL Test & ELISA Test	
Model	: VDS 180
Variable Speed	: 180 RPM

Specifications subject to change without notice

V. S. International Pvt. Ltd.

 A-204, Neelam Centre, Hind Cycle Road, Worli, Mumbai - 400 030., INDIA.
 Contact : +91 22 4333 3900 / 2498 8777 Fax: +91 22 6660 8779