

ACCUREX




TD-4230 Blood Glucose Monitoring System



STRIP
EJECTION



PC LINK



TD-4230

Blood Glucose Monitoring System

Product Features:

- Result indicator • Ketone warning
- 4 alarm settings • Internal code setting
- Alternative Site Testing (AST) • Strip-Ejection function
- 0.7 µl blood volume • 7 seconds reaction time
- 20 - 600 mg/dl measurement range
- 7, 14, 21, 28, 60 and 90 days average function
- 450 test memory with date and time
- PC download capability



Product Specifications:

Model Number	TD-4230
Result Range	20 - 600 mg/dl
Operating Condition	10°C - 40°C, below 85% R.H.
Sample Volume	0.7 µl
Reaction Time	7 seconds
Strip Ejection	Equipped
Ketone Warning	Enabled (≥ 240 mg/dl)
Alarm Setting	4 sets
Result Indicator	Enabled
Code Setting	Internal code setting 5/25 codes
Day Average	7, 14, 21, 28, 60 and 90 days average
Memory Capacity	450 (with date & time)
Power Source	1*3V CR2032
Device Size	88 mm (L) x 53 mm (W) x 15 mm (H)
Data Output	RS232/USB cable
Weight	43.58 gm.
Hematocrit Range	20 - 60%
Enzyme Type	GOD



ACCUREX BIOMEDICAL PVT. LTD.

Head Office - Mumbai. Tel. : 91(022) 6744 6744; Fax : 91(022) 6744 6755
E-mail : accurex@vsnl.com; Website : www.accurex.org
Plant : G-54, MIDC Tarapur, Boisar, Thane - 401506. India.

ACCUREX