

## All Age Digital Weighing Scale

(3 in 1)

200 Kg with +/- 20 gm Accuracy

- Durable, Moulded Plastic Baby Tray for Infant
- Folding Seat for Sitting Child
- Sturdy Platform for an Adult
- Separate LED Display Module to be kept on table or hanged on the wall
- Feather Touch Keys
- Zero Setting Facility
- Auto Hold Facility for Accurate Weighing of Active Babies
- Kg to Pound Conversion Facility
- Operating Voltage : 230V/50Hz AC



## All Age Digital Weighing Scale

(2 in 1)

200 Kg with +/- 20 gm Accuracy

- Baby Mattress for Infant
- Folding Acrylic Panel for Sitting Child
- Sturdy Platform for an Adult
- Feather Touch Keys
- Zero Setting Facility
- Auto Hold Facility for Accurate Weighing of Active Babies
- Kg to Pound Conversion Facility
- Operating Voltage : 230V/50Hz AC

